

Product Datasheet

Rat Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit (orb1664378)

Description	Rat Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 0.88 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.