

Product Datasheet

Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit (orb1664379)

Description	Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 0.96 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, urine or other biological fluids.
Expiration Date	Please enquire.