

## Product Datasheet

# Human Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit (orb1664380)

<b>Description</b>	Human Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 3.12-200 ng/mL. Sensitivity: 0.86 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, urine or other biological fluids.
<b>Expiration Date</b>	Please enquire.