

Product Datasheet

Rat Neurotensin (NT) ELISA Kit (orb1664385)

Description	Rat Neurotensin (NT) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4.68-300pmol/L. Sensitivity: 1.3pmol/L. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.