

## Product Datasheet

# Human N-Terminal Pro C-Type Natriuretic Peptide (NT-ProCNP) ELISA Kit (orb1664386)

<b>Description</b>	Human N-Terminal Pro C-Type Natriuretic Peptide (NT-ProCNP) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 15.6-1000 pg/mL. Sensitivity: 4.42 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
<b>Expiration Date</b>	Please enquire.