

Product Datasheet

PLD3 antibody (orb166442)

Description	Rabbit polyclonal antibody to PLD3 (N-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PLD3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54705
Uniprot ID	Q81V08
NCBI	NP_001026866.1 , NP_036400.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

