

Product Datasheet

Mouse Neuropeptide FF (NPFF) ELISA Kit (orb1664472)

Description	Mouse Neuropeptide FF (NPFF) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.048 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 μ L. Sample Type: urine, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.