

Product Datasheet

Canine Nitric Oxide Synthase 2, Inducible (NOS2) ELISA Kit (orb1664495)

Description	Canine Nitric Oxide Synthase 2, Inducible (NOS2) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.126 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.