

Product Datasheet

Human Neurokinin B (NKB) ELISA Kit (orb1664517)

Description	Human Neurokinin B (NKB) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5-800 pg/mL. Sensitivity: 3.12 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.