

Product Datasheet

Human Neurofilament 3 (NEF3) ELISA Kit (orb1664580)

Description	Human Neurofilament 3 (NEF3) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cerebrospinal fluid or other biological fluids.
Expiration Date	Please enquire.