

Product Datasheet

Porcine Myeloperoxidase (MPO) ELISA Kit (orb1664742)

Description	Porcine Myeloperoxidase (MPO) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.781-50 ng/mL. Sensitivity: 0.35 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.