

Product Datasheet

Human Multimerin 1 (MMRN1) ELISA Kit (orb1664758)

Description	Human Multimerin 1 (MMRN1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 9.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.