

Product Datasheet

Canine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit (orb1664844)

Description	Canine Macrophage Inflammatory Protein 1 Alpha (MIP1a) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 28 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.