

Product Datasheet

Human Maltase Glucoamylase, Intestinal (MGA) ELISA Kit (orb1664866)

Description	Human Maltase Glucoamylase, Intestinal (MGA) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.63 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.