

Product Datasheet

Human Myeloid Cell Leukemia Sequence 1, Bcl2 Related (MCL1) ELISA Kit (orb1664914)

Description	Human Myeloid Cell Leukemia Sequence 1, Bcl2 Related (MCL1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.118 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.