

Product Datasheet

Mouse Mannose Associated Serine Protease 1 (MASP1) ELISA Kit (orb1664935)

Description	Mouse Mannose Associated Serine Protease 1 (MASP1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2.34-150 ng/mL. Sensitivity: 0.8 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.