

Product Datasheet

Human Lymphotoxin Beta Receptor (LTbR) ELISA Kit (orb1664974)

Description	Human Lymphotoxin Beta Receptor (LTbR) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.104 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.