

Product Datasheet

Rat Low Density Lipoprotein Receptor Related Protein 8 (LRP8) ELISA Kit (orb1664990)

| | |
|----------------------------|--|
| Description | Rat Low Density Lipoprotein Receptor Related Protein 8 (LRP8) ELISA Kit |
| Reactivity | Rat |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 31.2-2000 pg/mL. Sensitivity: 13.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |