

Product Datasheet

Rat Lectin Like Oxidized Low Density Lipoprotein Receptor 1 (LOX1) ELISA Kit (orb1665027)

Description	Rat Lectin Like Oxidized Low Density Lipoprotein Receptor 1 (LOX1) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 pg/mL. Sensitivity: 22.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.