

Product Datasheet

CHORDC1 antibody (orb166511)

Description	Rabbit polyclonal antibody to CHORDC1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CHORDC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37490
Uniprot ID	Q9UHD1
NCBI	NP_036256.2 , NP_001137545.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

