

Product Datasheet

Mouse Integrin Alpha 3 (ITGa3) ELISA Kit (orb1665277)

Description Mouse Integrin Alpha 3 (ITGa3) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL. Sensitivity: 0.118 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates or other biological fluids.

Expiration Date Please enquire.