

Product Datasheet

Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit (orb1665285)

Description	Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 46.8-3000 pg/mL. Sensitivity: 16.13 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.