

Product Datasheet

Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit (orb1665286)

Description	Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.