

Product Datasheet

Bovine Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit (orb1665288)

Description	Bovine Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 46.8-3000 pg/mL. Sensitivity: 15.46 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.