

Product Datasheet

Rat Insulin Like Protein 3 (INSL3) ELISA Kit (orb1665333)

Description	Rat Insulin Like Protein 3 (INSL3) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.234-15 ng/mL. Sensitivity: 0.062 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.