

Product Datasheet

Mouse Interleukin 8 Receptor Beta (IL8Rb) ELISA Kit (orb1665367)

| | |
|----------------------------|---|
| Description | Mouse Interleukin 8 Receptor Beta (IL8Rb) ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids. |
| Expiration Date | Please enquire. |