

## Product Datasheet

# Mouse Interleukin 5 Receptor Alpha (IL5Ra) ELISA Kit (orb1665395)

<b>Description</b>	Mouse Interleukin 5 Receptor Alpha (IL5Ra) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 62.5-4000 pg/mL. Sensitivity: 28.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 $\mu$ L. Sample Type: serum, plasma or other biological fluids.
<b>Expiration Date</b>	Please enquire.