

Product Datasheet

Rat Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit (orb1665464)

Description	Rat Interleukin 22 Receptor Alpha 2 (IL22Ra2) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.105 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.