

Product Datasheet

Canine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit (orb1665489)

Description	Canine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 46.8-3000 pg/mL. Sensitivity: 14.35 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.