

Product Datasheet

Canine Interleukin 13 (IL13) ELISA Kit (orb1665575)

| | |
|----------------------------|--|
| Description | Canine Interleukin 13 (IL13) ELISA Kit |
| Reactivity | Canine |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Expiration Date | Please enquire. |