

Product Datasheet

Porcine Insulin Like Growth Factor 2 (IGF2) ELISA Kit (orb1665663)

Description	Porcine Insulin Like Growth Factor 2 (IGF2) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.64 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.