

Product Datasheet

Porcine Insulin Like Growth Factor 1 (IGF1) ELISA Kit (orb1665674)

Description	Porcine Insulin Like Growth Factor 1 (IGF1) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.146 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.