

Product Datasheet

Canine Insulin Like Growth Factor 1 (IGF1) ELISA Kit (orb1665679)

Description	Canine Insulin Like Growth Factor 1 (IGF1) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2.5-160 ng/mL. Sensitivity: 0.9 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.