

Product Datasheet

Porcine Immunoglobulin E (IgE) ELISA Kit (orb1665682)

Description Porcine Immunoglobulin E (IgE) ELISA Kit

Reactivity Porcine

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.066 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.