

Product Datasheet

Porcine Interferon Gamma (IFN γ) ELISA Kit (orb1665690)

| | |
|----------------------------|---|
| Description | Porcine Interferon Gamma (IFN γ) ELISA Kit |
| Reactivity | Porcine |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Expiration Date | Please enquire. |