

## Product Datasheet

# Guinea pig Interferon Gamma (IFN $\gamma$ ) ELISA Kit (orb1665691)

<b>Description</b>	Guinea pig Interferon Gamma (IFN $\gamma$ ) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 $\mu$ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.

---

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)