

## Product Datasheet

# Monkey Interferon Beta (IFN $\beta$ ) ELISA Kit (orb1665694)

<b>Description</b>	Monkey Interferon Beta (IFN $\beta$ ) ELISA Kit
<b>Reactivity</b>	Monkey
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.1 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 $\mu$ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.