

Product Datasheet

Human Interferon Gamma Inducible Protein 16 (IFI16) ELISA Kit (orb1665715)

Description	Human Interferon Gamma Inducible Protein 16 (IFI16) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.049 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.