

Product Datasheet

Rat Insulin Degrading Enzyme (IDE) ELISA Kit (orb1665726)

Description	Rat Insulin Degrading Enzyme (IDE) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.129 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.