

Product Datasheet

Mouse Insulin Degrading Enzyme (IDE) ELISA Kit (orb1665727)

Description	Mouse Insulin Degrading Enzyme (IDE) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 pg/mL. Sensitivity: 1.33 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.