

Product Datasheet

Human 5-Hydroxytryptamine Receptor 4 (HTR4) ELISA Kit (orb1665750)

Description	Human 5-Hydroxytryptamine Receptor 4 (HTR4) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.063 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.