

Product Datasheet

Mouse Heat Shock 70kDa Protein 8 (HSPA8) ELISA Kit (orb1665785)

Description	Mouse Heat Shock 70kDa Protein 8 (HSPA8) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.116 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Expiration Date	Please enquire.