

Product Datasheet

Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit (orb1665794)

Description	Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.14 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.