

Product Datasheet

Human Heat Shock Protein 90kDa Alpha A1 (HSP90aA1) ELISA Kit (orb1665810)

| | |
|----------------------------|--|
| Description | Human Heat Shock Protein 90kDa Alpha A1 (HSP90aA1) ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 1.56-100 ng/mL. Sensitivity: 0.55 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Expiration Date | Please enquire. |