

Product Datasheet

Mouse Haptoglobin (Hpt) ELISA Kit (orb1665846)

Description	Mouse Haptoglobin (Hpt) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 ng/mL. Sensitivity: 13.9 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.