

Product Datasheet

Canine Haptoglobin (Hpt) ELISA Kit (orb1665847)

Description	Canine Haptoglobin (Hpt) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 ng/mL. Sensitivity: 24.5 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.