

Product Datasheet

Mouse High Mobility Group Protein 1 (HMG1) ELISA Kit (orb1665872)

| | |
|----------------------------|--|
| Description | Mouse High Mobility Group Protein 1 (HMG1) ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 62.5-4000 pg/mL. Sensitivity: 23.73 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids. |
| Expiration Date | Please enquire. |