

Product Datasheet

Canine High Mobility Group Protein 1 (HMG1) ELISA Kit (orb1665874)

Description	Canine High Mobility Group Protein 1 (HMG1) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.24 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.