

Product Datasheet

Human Hemoglobin Gamma 1 (HBg1) ELISA Kit (orb1665935)

Description	Human Hemoglobin Gamma 1 (HBg1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 468-30000 ng/mL. Sensitivity: 122.6 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or erythrocyte lysates.
Expiration Date	Please enquire.