

## Product Datasheet

# Human Hemoglobin Beta (HBb) ELISA Kit (orb1665940)

**Description** Human Hemoglobin Beta (HBb) ELISA Kit

**Reactivity** Human

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 1.56-100 µg/mL. Sensitivity: 0.45 µg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or erythrocyte lysates.

**Expiration Date** Please enquire.